


## CO<sub>2</sub> (capnophilic) cultivation

For the cultivation at approximately 5% of CO<sub>2</sub>

AnaeroPack®-CO<sub>2</sub>  
AnaeroPouch®-CO<sub>2</sub>

Will support the growth of *Hemophilus* spp. and *Neisseria* species.

### Atmospheric profile by AnaeroPack®-CO<sub>2</sub>



AnaeroPack®-CO<sub>2</sub>



AnaeroPouch®-CO<sub>2</sub>

	Rect. Jar			W-Zip Pouch or Pouch-Bag
	2.5L	7L	0.4L	
AnaeroPack®-CO <sub>2</sub>	1 sachet	3 sachets		
AnaeroPouch®-CO <sub>2</sub>			1 sachet	1 sachet ✳

✳ Please put 2 plates. For the culture of only one plate, please add one non-inoculated plate to maintain the proper environment inside the bag.