

“

INNOVITA is a leading manufacturer of diagnostic solutions for healthcare, striving for a more efficient healthcare system to enhance the health and well-being of people all over the world.

”

INNOVITA BIOLOGICAL TECHNOLOGY CO., LTD.

Head Office (Beijing):

Block 6, Courtyard No. 1, Haiying Road, Fengtai District, Beijing, China 100070


Branch (Tangshan):

No. 699, Juxin Street, High-Tech Industrial Development Zone, Qian'an, Hebei, China 064400

Tel: +86 10 83681277

Fax: +86 10 83682966

Email: export@innovita.com.cn



Daya Suar Scientific

2019-nCoV

Ag Test (Nasopharyngeal)

For Professional Use



Innovita 2019-nCoV Antigentest is intended for the direct and qualitative detection of SARS-CoV-2 nucleocapsid protein antigen in nasopharyngeal swabs for screening of individuals without symptoms or other reasons to suspect COVID-19 infection. Easily test for SARS-CoV-2 infection whether you're symptomatic or asymptomatic.

Features



Convenient

Suitable for testing in areas where laboratory testing is not available (Professional use)



Time-saving

Getting quick result within 15-30 minutes



Easy-handling

No special training or device



Reliable

Better clinical performance



+62 21 38777055
+62 811 811 87 857
sales@dayasuarscientific.co.id
www.dayasuarscientific.co.id

PT Daya Suar Scientific

Rukan Botanic Mega Kebon Jeruk Blok H.8 No.7,
Desa/Kelurahan Joglo, Kec. Kembangan,
Jakarta Barat 11640


Partnering For Excellence

How It Works

● Specimen Collection



● Specimen Handling



Result Interpretation



Positive



Positive



Negative



Invalid



Invalid



Clinical performance

Sensitivity

98.7%

Specificity

100%

INNOVITA 2019-nCoV Ag Test	Number of PCR Positive samples	
	Positive	Negative
Positive	148	0
Negative	2	250
total	150	250

Relative Sensitivity: 98.7% (95% CI: 95.27% to 99.84%)

Relative Specificity: 100% (95% CI: 98.53% to 100%)

*Confidence interval

Specifications

Packing Size: 1/25 tests /box

Sample Type: nasopharyngeal swab specimen

Storage: 4°C -30°C

Shelf Life: 18months